





MFM imaging with A-Probe

Topography



Magnetic



Sample : Hard Disk
Scan size: 15 μm x 15 μm
Lifted height : 40 nm

Akiyama-Probe with magnetic coating*

on Veeco Nanoscope III, with Nanosurf easyPLL, at IMT Univ. of Neuchâtel Switzerland

*This short note is to demonstrate a potential of A-Probe. A-Probe with magnetic coating is not yet optimized and not ready for distribution.

NANOSENSORS™
Neuchatel Switzerland
<http://www.nanosensors.com>
e-mail: info@nanosensors.com